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$$\text{Free carnitine concentration (}\mu\text{mol/L)} = \frac{\left(\text{Absorbance of the sample} \right) - \left(\text{Absorbance of the reagent blank} \right)}{\left(\text{Absorbance of the standard} \right) - \left(\text{Absorbance of the reagent blank} \right)} \times \text{Concentration of the standard (}\mu\text{mol/L)}$$

Free Carnitine "Kainos"

This product is warranted to perform as described in its labeling and literature when used in accordance with all instructions.

【Performance】

- Sensitivity and measurement range
 - Change in the OD of saline as the sample is within 0.02/minute.
 - Change in the OD of 50 $\mu\text{mol/L}$ standard solution as the sample is in the range of 0.02-0.03/minute.
 - The measurement range is 0.5-200 μmol of carnitine/L.
- Specificity
The specificity is within 90-110% of the known concentration sample.
- Reproducibility
The CV of measured values is less than 5.0%.
- Correlation
This Free Carnitine "Kainos" was evaluated against Kainos's previous free carnitine measurement kit. With 50 serum samples, a good correlation with $r=0.997$ and the good regression formula with $y = 0.980x - 0.652$ were obtained.

【HANDLING and OPERATION PRECAUTION】

- If serum sample is used, use fresh serum even though other substances contained in serum have little effect on the assay.
- Anticoagulants that are within normal usage concentration ranges do not affect the assay.
- Preservatives are within normal usage concentration ranges do not affect the assay except for 5% formalin.
- Make sure to use reagents that are from the same lot.
- Before performing an assay, rinse the measurement apparatus well and calibrate using both blank and standard.
- In case carnitine concentration in the sample exceeds the measurable range, dilute the sample with saline and repeat the assay.
- Handling of human blood:
Serum or other blood samples should be handled as potentially infectious as they may contain infectious agents including HB antigen. Any instruments that have come to contact samples should also be treated potentially infectious.
- Handling of reagents that contain sodium azide:
Reagents in this kit contain sodium azide as preservatives. Because sodium azide may react with metals including lead and copper to produce explosive metal azides, any waste solution should be flushed with copious amount of water. In case any of the reagents contact eyes, mouth and/or skin, immediately flush with copious amount of water then consult a physician as necessary.
- This kit contains two reaction mixtures. Care must be taken during measurement not to mix up these reagents.
- The kit should be used before the expiration date printed on the product label.
- The kit is intended only for research purposes.

【WARRANTY】

This product is warranted to perform as described in its labeling and literature when used in accordance with all instructions. KAINOS LABORATORIES, Inc. DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, and in no event shall KAINOS LABORATORIES, Inc. be liable for consequential. Replacement of the product or refund of the purchase price is the exclusive remedy for the purchaser. This warranty gives you specific legal right and you may have other rights which vary from state to state.

【REFERENCES AND CUSTOMER SERVICE】

- REFERENCES
 - Takahashi, M., et al: Clin.Chem., 40: 817 (1994)
 - Tsuji, A: Saibokogaku, 18: 11 (1999) (in Japanese)
 - Kawano, N: Nihon Rinsho, 57: Special Edition (1999) (in Japanese)

2. CUSTOMER SERVICE

Kainos Laboratories, Inc.
38-18, Hongo 2-Chome Bunkyo-ku, Tokyo, Japan 113-0033

Manufacturer

K KAINOS LABORATORIES, INC.

38-18 Hongo 2-Chome Bunkyo-ku, Tokyo 113-0033 Japan

In technical cooperation

AsahiKASEI
ASAHI KASEI PHARMA

1-105 Kanda Jinbocho Chiyoda-ku, Tokyo 101-8101 Japan

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
It is well known fatty acids are utilized as energy sources by being broken down through the β -oxidation cycle to generate ATP. Since this reaction takes place within mitochondria, fatty acids must be transported into mitochondria. The transportation of fatty acid may not be direct because the mitochondrial membranes are not permeable to free fatty acids. Fatty acids require some carrier to enter mitochondria, which is the role of carnitine in the mitochondrial membranes.

The Kainos-developed Free Carnitine "Kainos" is a reagent with superior stability and ease of use that measures the level of free carnitine based on an enzymatic cycling method.

【PRINCIPLE OF THE ASSAY】

The assay principle is enzymatic cycling method.

Carnitine in a given sample is specifically oxidized by Carnitine Dehydrogenase (CDH) under the presence of oxidized β -Thionicotinamide Adenine Dinucleotide (Thio-NAD⁺) to dehydrocarnitine and reduced β -Thionicotinamide Adenine Dinucleotide (Thio-NADH). Dehydrocarnitine in turn is specifically reduced by CDH under the presence of reduced β -Nicotinamide Adenine Dinucleotide (NADH) to carnitine and oxidized β -Nicotinamide Adenine Dinucleotide (NAD⁺). The following enzymatic cycling reaction leads to the accumulation of Thio-NADH. Because the generated amount of Thio-NADH is proportional to the carnitine concentration in the sample, free carnitine concentration can be determined by the specific absorption of Thio-NADH.

【FEATURE】

- High sensitivity and accuracy based on an enzymatic cycling method.
- High substrate specificity of the enzyme used.
(Unwanted interference by other similar substances is eliminated.)
- Adaptable to automated systems.
(Large number of samples can be quickly processed.)

【REAGENTS (Catalog No. RK-2400)】

- | | |
|---|--------------------|
| 1) Reagent 1 Thio-NAD ⁺ , Sodium Azide | :1 bottle(20 mL) |
| 2) Reagent 2A Sodium Azide | :1 bottle(10 mL) |
| 3) Reagent 2B Carnitine Dehydrogenase (CDH) and NADH | :1 vial(for 10 mL) |
| 4) Standard Solution 50 $\mu\text{mol/L}$ L-Carnitine and Sodium Azide | :1 bottle(5 mL) |

* These reagents should be stored at 2-10 °C in dark.

【PREPARATION of REAGENT】

- Reagent 1(R-1)
No preparation is required.
Keep lid closed when not in use. May be kept for 1 month at 2-10 °C in dark.
- Reagent 2(R-2)
Dissolve reagent 2B in reagent 2A.
Keep lid closed when not in use. May be kept for 1 month at 2-10 °C in dark.

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- Measurement (Example: Roche HITACHI 917)